# IKA® EH 4 basic

Distributed by:

**Lab Unlimited**  
CARL STUART GROUP

Tallaght Business Park  
Whitestown, Dublin 24,  
Ireland  
D24 RFK3

Tel: (01) 4523432  
Fax: (01) 4523967  
E-mail: [info@labunlimited.com](mailto:info@labunlimited.com)  
Web: [www.labunlimited.com](http://www.labunlimited.com)

Quatro House, Frimley Road,  
Camberley,  
United Kingdom  
GU16 7ER

Tel: 08452 30 40 30  
Fax: 08452 30 50 30  
E-mail: [info@labunlimited.co.uk](mailto:info@labunlimited.co.uk)  
Web: [www.labunlimited.co.uk](http://www.labunlimited.co.uk)

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## CE-DECLARATION OF CONFORMITY

EN

We declare under our sole responsibility that this product corresponds to the regulations 89/336/EEC and 73/23/EEC and conforms with the standards or standardized documents DIN EN IEC 61010-1 and DIN EN IEC 61326-1

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## Garantie

You have purchased a laboratory machine which meets the highest engineering and quality standards.

In accordance with guarantee conditions, the guarantee period is 24 months. For claims under the guarantee please contact your local dealer. You may also send the machine direct to our works, enclosing the delivery invoice and giving reasons for the claim. You will be liable for freight costs.

The warranty does not cover wearing parts, nor does it apply to faults resulting from improper use or insufficient care and maintenance contrary to the instructions in this operating manual.

If you require servicing, return the appliance in its original packaging. Storage packaging is not sufficient. Please also use suitable transport packaging.

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# Contents

	Page
Guarantee	2
Safety instructions	8
Proper use	9
Unpacking	10
Commissioning	10
Installation, filling and emptying	10
Starting	10
Setting of the safety temperature	11
Maintenance and cleaning	11
Associated standards and regulations	11
Accessories	12
Technical data	12
List of spare parts EH4 basic	12
Spare parts diagram EH4 basic	19

# Safety instructions


**Caution:** With this equipment only media may be worked on/or heated, whose flash point lies over the adjusted safety temperature delimitation. Value for safety temperature limit see chapter „Technical data“. Examine regularly the function of the safety temperature limiter (see section “Setting of the safety temperature”).

The instrument may not be operated in rooms with explosion hazard.

Never place the unit in operation without sufficient bath liquid.

**Risk of burns:** Danger when touching the housing parts and the heater coils. The temperature of the heater coils can very rise. Before placing the device in service, please make certain that it is set to the lowest tem-

perature. Otherwise, the device will begin heating up to the last temperature to be set.

After an interruption in the power supply during the temperature control process, the device starts back up by itself.

Coverings or parts that can be removed from the device without the aid of any tools must be put back on the device again to ensure safe operation, for example to keep foreign objects and liquids, etc. from getting into the device.

Please observe applicable safety instructions and regulations as well as work protection and accident prevention requirements for use in the laboratory.

The device must be fastened onto the laboratory stand so that it cannot tip over. The support brackets on the bath vessel (included with delivery) or the boom, which screws into the threaded socket, can be used for this purpose. The device is not designed for manual operation.

You should also be careful to ensure that the minimum clearances and immersion depths in the temperature control medium are observed (*see diagram*).

Exercise caution while filling up hot bath vessels.

Do not transport or empty the bath vessel while it is still hot. This may result in accidents, especially scalding injuries.

If you are using plastic bath vessels, observe the limited range of working temperatures and use a suitable bath liquid (it may be necessary to consult with the bath manufacturer).

The safety circuit (safe temperature) must be adjusted so that the maximum permissible temperature cannot be exceeded even in the event of a fault (see the section "Adjusting the safety temperature").

When working at elevated temperatures, parts of the bath vessel may take on high surface temperatures. Exercise caution when coming in contact with these surfaces.

Secure distributor hoses and tubes against slipping and avoid kinks. Check distributor hoses, distributor tubes and bath vessels at regular intervals for possible material fatigue (cracks/leaks).

Demineralized water is recommended as a bath or temperature control liquid.

IKA assumes no liability if temperature control liquids not recommended by IKA are used.

**The maximum temperature for water as a bath liquid/temperature control liquid is 90°C.** If water is used as a bath liquid above this temperature, there is heavy loss of liquid due to evaporation, associated with steam condensation in and on parts in the area of evaporation.

There is a risk of damage being caused by the bath running dry quickly. Related to this, damage may also be caused to the EH4 basic by water condensation.

Loss of the bath liquid resulting from evaporation and condensation of water vapor on the EH4 basic can be kept to a minimum by using suitable bath coverings.

Continuous monitoring of the bath and the filling level of the bath liquid is required, especially at high temperatures.

The device must only be operated with grounded mains power connections.

Please be aware of the possible danger posed by toxic (poisonous) or flammable gasses being liberated from material that is being heated up.

## Proper use

The EH4 basic suspended thermostat is suitable for use in modifying the temperature of liquids from room temperature to 100 °C in open containers.

## Unpacking

Please unpack the equipment carefully and check for any damages. It is important that any damages which may have arisen during transport are ascertained when unpacking. If applicable a fact report must be set immediately (post, rail or forwarder).

The delivery scope covers: an EH4 basic and operating instructions.

## Commissioning

Check whether the voltage specified on the type plate matches the mains voltage available. The power socket used must be earthed (protective earth conductor contact). If these conditions are met, the device is ready to operate after plugging in the mains plug. If these procedures are not followed, safe operation cannot be guaranteed and/or the equipment may be damaged. Observe the ambient conditions (temperature, humidity, etc.) listed under Technical Data.


## Installation, filling and emptying

Place the IKA EH4 basic suspended thermostat on the bath vessel and secure the EH4 basic suspended thermometer with the fixing clamp on the bath vessel.

Fill the bath vessel with liquid to moderate the temperature (when recommend water with moderate mineral content). Make certain the bath is filled to the level indicated.

To empty the bath vessel, you should remove the ET basic suspended thermostat and all accessory pieces from the bath vessel.

## Starting

The EH4 basic suspended thermostat is turned on with the rocker switch **(A)**. The „Temp“ rheostat dial **(B)** with accompanying temperature of the bath.

The heating-up phase is indicated by the yellow pilot light **(C)**. When the required temperature is reached, it flashes on and off with the rhythm of the heater.

The bath temperature that was set can be verified with the monitoring thermometer.

## Setting of the safety temperature

The desired safety temperature can be adjusted with the „safe temp“ set screw. As soon as the desired safety temperature is exceeded by a higher target temperature that has been set, or if an error is encountered in the temperature controller, the red pilot light **(E)** „Alarm“ lights up. The heater turns off permanently in this state, and the pump continues to run.

The unit can only be placed in operation again by switching off and then back on again with the rocker switch **(A)**.

After the target temperature has been reached, the safety temperature is determined as follows:

- Using an appropriate tool, turn the „safe Temp“ set screw **(D)** of the safety temperature controller until the switch-off point to the left. The red pilot light **(E)** „Alarm“ lights up.
- Then turn the „safe Temp“ set screw **(D)** of the safety temperature controller just slightly to the right.
- After you have turned the unit off and back on again with the rocker switch, it is ready for operation.

## Maintenance and cleaning

The IKA EH4 basic is maintenance-free. It is subject only to the natural wear and tear of components and their statistical failure rate.

**When ordering spare parts, please give the manufacturing number shown on the type plate, the machine type and the name of the spare part.**

Please send in equipment for repair only after it has been cleaned and is free from any materials which may constitute a health hazard. Use only cleansing agents which have been approved by IKA to clean IKA devices. To remove use:

Dyes	isopropyl alcohol
Construction materials	water containing tenside / isopropyl alcohol
Cosmetics	water containing tenside / isopropyl alcohol
Foodstuffs	water containing tenside
Fuels	water containing tenside

For materials which are not listed, please request information from IKA. Wear the proper protective gloves during cleaning of the devices.

Electrical devices may not be placed in the cleansing agent for the purpose of cleaning.

Before using another than the recommended method for cleaning or decontamination, the user must ascertain with the manufacturer that this method does not destroy the instrument.

To avoid contamination, it is essential to change the bath liquid from time to time.

## Associated standards and regulations

**Construction in accordance with the following safety standards**

DIN EN IEC 61010-1 UL 3101-1 CAN/CSA C22.2 (1010-1)

**Construction in accordance with the following EMC standards**

DIN EN IEC 61326-1

**Associated EU guidelines**

EMC-guidelines: 89/336/EG

Machine guidelines: 73/23/EG

## Accessories

<b>EH 4.1</b>	Bath vessel
<b>EH 4.2</b>	Bath vessel
<b>EH 4.3</b>	Bath vessel

## Technical data (acc. DIN 12876)

design voltage:	<b>VAC</b> 230±10%
or	<b>VAC</b> 115±10%
design frequency:	<b>Hz</b> 50/60
class partitioning	I
acc. DIN 12876	
identifying marking	NFL
acc. DIN 12876	
heating power:	<b>W</b> 1500
working temperature range:	<b>°C</b> 25 to 100
temperature setting:	infinitely variable from ambient temperatur
temperature constancy:	<b>K</b> ±0,12 by 70°C and 12 Ltr H <sub>2</sub> O
temperature control:	two-step control
control sensor:	electro-mechanical
safety sensor:	electro-mechanical
safety temperature- limitation adjustable:	<b>°C</b> 25 to 200
pump type:	pressure pump
max. pump pressure:	<b>bar</b> 0,08 (0 litre flow of discharge)
max. flow of discharge:	<b>l/min</b> 5 (0 bar back pressure)
min. bath depth:	<b>mm</b> 160
min. effective depth:	<b>mm</b> 75
permitted ambient temperature:	<b>°C</b> +5 to +40
permitted relative humidity:	<b>%</b> 80

permitted duration of operation:	<b>%</b> 100
protection class acc. DIN 60529:	IP 31
Fuses:	230V <b>A</b> 10AT Id.Nr. 27 554 00
	115V <b>A</b> 15AT Id.Nr. 26 494 00
contamination level:	2
overvoltage category:	II
protection class:	I (protective earth)
Operation at a terrestrial altitude:	<b>m</b> max. 2000 above sea level
dimensions: (W x D x H)	<b>mm</b> 105 x 139 x 319
weight:	<b>kg</b> 2,3

## List of spare parts EH4 basic

Pos.	Designation	Pos.	Designation
1	Plastite-Counter sunk screw	36	Spacer piece
3	Plastite-Cylinder head screw	37	Suspension
6	Hexagonal nut	38	Housing-lower part
8	Disk	39	Pump housing
9	Snap ring	40	Pump blade
11	Connection line	41	Pump cover
12	Strain relief	42	Holding angle
13	Pilote lamp yellow	43	BLP - power supply
14	Pilote lamp red	44	Housing-upper part
15	Rocker switch	45	Front panel
16	Clamp	46	Operatin knob
17	Clamp		
20	Safety capillary tube controller		
21	Capillary tube controller		
28	Motorholder	47	Cap
29	Bottom plate	54	Covering
30	Insulating disk	60	Container clamp
31	Pump holder	2001	Shaded pole motor,
34	Protective tube		compl.
35	Tube heater body	5001	Fuse

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